
BioSerenity
Medical Devices Group

Wearable Products Line

**USER MANUAL / MANUEL UTILISATEUR /
MANUALE UTENTE / GEBRAUCHSANWEISUNG**


Table of Contents

1. FOREWORD	6
1.1. Intended Use	6
1.2. Indications for Use	6
1.3. Scope	6
1.4. Intended Benefits	6
1.5. Document History	7
1.6. Glossary	7
1.7. Abbreviations	7
1.8. Symbology	7
2. INTRODUCTION	8
2.1. Contraindications	8
2.2. Adverse Reaction (for IceCap Neonate models)	8
2.3. System Overview	8
2.3.1. IceCap 2 and IceCap 2 Small	10
2.3.2. IceCap Neonate	10
2.3.3. Neuronaute® Adapters	10
2.4. Compatible Devices	11
3. APPLYING THE ICECAPS	12
3.1. Checking the Device Before Application	12
3.1.1. Checking the Device	12
3.1.2. Checking the patient	12
3.1.3. Checking the Environment	13
3.1.4. Size Chart	13
3.2. Applying the IceCap 2	13
3.2.1. Required Equipment	14
3.2.2. Applying the IceCap 2	15
3.3. Applying the IceCap Neonate	18
3.3.1. Required Equipment	18
3.3.2. Applying the IceCap Neonate	19
4. SAFETY INFORMATION	23
4.1. Operating Conditions	23
4.2. Safety Information	23
4.3. Information on Residual Risks	24
5. REGULATORY INFORMATION	26
5.1. System Life Cycle	26
5.1.1. Product Lifetime	26
5.1.2. Recycling the IceCap 2 Parts	26
5.2. Product Symbology	27

5.2.1. Manufacturer Information 28

5.3. Medical Device Classification 28

5.4. Contact 28

5.4.1. Contact Information 28

5.4.2. Internet 28

5.4.3. Instructions for Use 28

5.4.4. Materiovigilance 28

1. FOREWORD

The **Wearable Products Line (WPL)** is a range of medical devices.

1.1. Intended Use

The **IceCaps** are medical devices used as EEG electrodes. They are used by Healthcare Professionals on a patient in case of neurological disorders with a short-term or long-term EEG record (up to 72 hours).

IceCap 2 shall be placed on patients weighing at least 10 kg (22.05 lbs) and having a head circumference above 43 cm (16.93 inches).

IceCap Neonate shall be placed on the head of babies, newborns and premature babies.

1.2. Indications for Use

Products of the **WPL** are intended to be used by a healthcare professional on a patient.

1.3. Scope

The purpose of this manual is to explain the features of the **WPL** to the user.

The **WPL** is composed of the following elements:

- IceCap 2
- IceCap 2 Small
- IceCap Neonate

The **IceCap Neonate** range is composed of a device that exists in three different sizes:

- IceCap Neonate (XS)
- IceCap Neonate (S)
- IceCap Neonate (M)

In this document, the term **IceCap Neonate** includes all products from the range.

1.4. Intended Benefits

The benefits of the device are:

- Collection of high quality EEG signals.
- Better use for healthcare professionals:
 - Electrodes are preplaced at the right position, in accordance with the international 10-20 system.
 - Reduction of the setup time on a patient.

1.5. Document History

The following table shows the document changes history since its creation.

Version	Date	Description
A	February 23 2023	Initial creation
H	March 18 2025	Updates and modifications

1.6. Glossary

Lead

An electrode used in an electroencephalogram (EEG) which contributes to a graphic representation of electrical brain activity.

1.7. Abbreviations

Abbreviation	Term
EEG	ElectroEncephaloGram
EMC	ElectroMagnetic Compatibility
ICE	Instant Care EEG
WEEE	Waste Electrical and Electronic Equipment
WPL	Wearable Products Line

1.8. Symbology

Note

Gives additional information for optimal use of the product.

Warning

Draws attention to important information and / or action: damage to the equipment is possible.

2. INTRODUCTION

2.1. Contraindications

Do not use the system in any of the following cases:

- Sensory processing disorder having an effect on skin sensitivity.
- Behavioral problems making the patient excessively agitated or aggressive.
- Mental disorder that is incompatible with using the device.
- Open wound in an area either covered or surrounded by the device.
- Use of the data gathered by the device nor any of its components to establish a brain death diagnosis.
- Allergy to one of the components of the device: silver, silver chloride, silicone, acrylic, polyethylene terephthalate (PET), polyethylene terephthalate glycol (PETG), polyether polyurethane, oxirane, (phenoxymethyl), acrylic oligomer, isobornyl acrylate esters, 2,2-dimethoxy-2-phenylacetophenone, benzophenone, talc, any synthetic material.
- The **IceCap 2** and **IceCap Neonate** are not to be used for electrical stimulation.

2.2. Adverse Reaction (for IceCap Neonate models)

Due to the nature of newly born infant skin, there may be some adverse reaction at the electrode sites. Thus, it is recommended that the skin at or around the electrode contact sites be monitored for any signs of skin reaction.

2.3. System Overview

The system is composed of the pieces of the equipment listed below.

Warning

The use of any unit not recommended in this document could have an unwanted effect on the system.

Note

Figures from this document may vary slightly from the actual device depending on the version used.

IceCap 2

It enables up to 72 hours of recording.


It consists of:

- An IceCap 2
OR
An IceCap 2 small
- Consumables (not provided)

It is used along with:

- An IceAdapter
OR
A Touchproof Adapter
- A DB25 cable

For more information, refer to [Section 3.2, “Applying the IceCap 2”](#).

IceCap Neonate

The **IceCap Neonate** enables up to 72 hours of recording.


It consists of:

- An IceCap Neonate
- Consumables (not provided)


It is used along with:

- An IceAdapter
OR
- A Touchproof Adapter
- A DB25 cable

For more information, refer to [Section 3.3, “Applying the IceCap Neonate”](#).


2.3.1. IceCap 2 and IceCap 2 Small


The **IceCap 2** is a single use cap which connects to a compatible EEG acquisition system using an IceAdapter or a Touchproof Adapter (refer to [Section 3.2.2, “Applying the IceCap 2”](#)).


2.3.2. IceCap Neonate

XS

S

M



The **IceCap Neonate** is a single use cap which connects to a compatible system using an IceAdapter or a Touchproof Adapter (refer to [Section 3.3.2, “Applying the IceCap Neonate”](#)).

In this document, the term **IceCap Neonate** includes all products from the range.

2.3.3. Neuronaute® Adapters



The Adapters are accessories that enable connecting a cap to a compatible EEG acquisition device using a specific cable.

IceAdapter

The IceAdapter connects to a compatible device using a DB25 cable.

Touchproof Adapter

The Touchproof Adapter connects to a compatible EEG acquisition device using a DB25 into 21 touchproof connectors cable.

2.4. Compatible Devices

The **WPL** is intended to be used with:

- **The Neuronaute system**

For more information regarding the **Neuronaute®** system, refer to Bioserenity documentation.

- Micromed recorder equipped with appropriate female DB25 connector using **IceAdapter**.
- Recorders equipped with female touchproof connector DIN 42802 using **TouchProof Adapter**.

3. APPLYING THE ICECAPS

3.1. Checking the Device Before Application

3.1.1. Checking the Device

Before using the system, check the device:

- Check that there is no damage on any part of the system.
- Check that the system is not damp, wet or dirty.

Warning

The IceCaps are a single use devices. They shall not be washed.

Warning

Never use a damaged or degraded device. In case of damaged or degraded device, contact the support.

Warning

Avoid tightly folding the tracks or connector during application.

3.1.2. Checking the patient

Before using the system, verify the patient condition and make sure that the patient has no contraindications.

The application shall be done on still patients (example: once the patient is asleep).

Warning

If the system is painful to wear, remove the device and contact a healthcare professional.

Warning

If skin redness occurs, remove the device and contact a healthcare professional.

3.1.3. Checking the Environment

Before using the system, note that the user should:

- Keep away from fire or any heat source that may cause damage to the device.
- Keep the device away from any source of water while recording.
- Remove the device in case of excessive sweating.
- Do not smoke near to the device.

Warning

Make sure that the environment is dry and there is no risk of the system getting wet while recording. In the event of excessive sweating, it is recommended to remove the system.

3.1.4. Size Chart

The **IceCaps** exist in different sizes. Before using the device, make sure to determine the correct size using the size chart below.

To measure the head perimeter of the patient, place a tape measure on the protuberance located in the back of the head and on the forehead.

If the patient is between two sizes, choose the inferior one.

	Minimum head perimeter [cm/(inch)]	Maximum head perimeter [cm/(inch)]
IceCap Neonate (XS)	23 (9.06)	31 (12.21)
IceCap Neonate (S)	30 (11.82)	38 (14.97)
IceCap Neonate (M)	37 (14.57)	51 (20.08)
IceCap 2 Small	43 (16.93)	53 (20.87)
IceCap 2	53 (20.87)	60 (23.63)

3.2. Applying the IceCap 2

The IceCap 2 records:

- 19 EEG channels (Fp1, Fp2, F7, F3, Fz, F4, F8, C3, Cz, C4, P3, Pz, P4, T3, T4, T5, T6, O1, O2)
- The Oz electrode is used for ground connection.
- The FpZ electrode is used as a reference for EEG calculation.

Warning

The IceCap 2 is a single use device.

The IceCap 2 cleaning and re-usability has not been validated and could impact the patient safety and the performances of the device.

Warning

The position of electrodes (IceCap 2 and IceCap 2 Small) can shift during the wear duration. Please ensure that the device is properly donned and positioned for the duration of use. External protection such as soft bandages or wraps should be used to maintain the positioning of the device for the duration of use. Failure to do so can result in signal capture deterioration.


3.2.1. Required Equipment

Before applying the **IceCap 2**, make sure that the following components are available.


Warning

The conductive paste must be designed for EEG electrodes.


Conductive paste and other creams must be used under the operating conditions provided by the manufacturer.

IceCap 2


Gloves


Swabs


Compatible system

Skin preparation gel
(mandatory if recording
duration > 12h)


Conductive paste


Adhesive cream

Cups x2

Medical tape and hair net

Pads

3.2.2. Applying the IceCap 2

Preparing the Application

Before applying the **IceCap 2**, make sure that all the components described previously are available and that the patient has no contraindications.


To prepare the application:

- Put some skin preparation gel in a cup.
- Put some conductive paste in a cup.
- Put the gloves on.


Applying the IceCap 2

To apply the **IceCap 2**:

1. Stretch the IceCap 2 and put it on the head of the patient:
 - Peel the protective film off from the **Front Nz** strip and fix it to the forehead (between the eyebrows), using the Cz electrode as a visual reference for the midline central location of the head.


- Peel the protective film off from the **Left** and **Right** strips and fix them to the temples.

2. Peel the protective film off from the **Back** strip and fix it to the back of the neck, below the hair line.


**Warning**

Fix the connectors using medical tape or a hair net to avoid any risk of strangulation

**Warning**

Make sure that the conductive side is in contact with patient's skin. Labels shall be on the outer side.

3. Check the electrodes positioning.
4. Using the wooden stick, remove the hair under each electrode and apply some skin preparation gel to the scalp.
5. Using the swabs, apply some conductive paste under each electrode.
6. Plug the **IceCap 2** to the IceAdapter or Touchproof Adapter.

**Warning**

Make sure that the L and R marks from the flat cable correspond to the figure above.

**Note**

Other compatible devices are listed in [Section 2.4, “Compatible Devices”](#).

7. Plug the IceAdapter or Touchproof Adapter to a compatible device.
8. Using a compatible EEG acquisition device, check the impedance. If the values are too high, remove the hair under the corresponding electrodes and softly push the electrodes against the scalp.
9. If additional fixation for electrode is needed, apply some adhesive cream on each pad square and place each pad square on top of electrode to help fixing the electrode on the skin.
10. Place each pad square on an electrode, with the adhesive cream in contact with the electrode. Softly push the pad to spread the cream.
11. Check the impedance as described in step 8.

3.3. Applying the IceCap Neonate

The **IceCap Neonate** is composed of:

- 11 EEG channels (size M)
or
9 EEG channels (size XS, S)
- 1 ground connection
- 1 reference electrode for EEG calculation.

Warning

The **IceCap Neonate is a single use device**

The **IceCap Neonate cleaning and re-usability has not been validated and could impact the patient safety and the performances of the device.**

Warning

The position of electrodes (IceCap Neonate**) can shift during the wear duration. Please ensure that the device is properly donned and positioned for the duration of use. External protection such as soft bandages or wraps should be used to maintain the positioning of the device for the duration of use. Failure to do so can result in signal capture deterioration.**

Warning

The patient must be under supervision during the record.

3.3.1. Required Equipment


Before applying the **IceCap Neonate**, make sure that the following components are available.

Warning


The conductive paste must be designed for EEG electrodes.


Conductive paste and other creams must be used under the operating conditions provided by the manufacturer.

IceCap Neonate



Gloves


Medical tape and hair net


Compatible system

Skin preparation gel
(mandatory if recording
duration > 12h)


Conductive paste

Adhesive cream

Pads

Cups x2

Swabs

3.3.2. Applying the IceCap Neonate

Preparing the Application

Before applying the **IceCap Neonate**, make sure that all the components described previously are available and that the patient has no contraindications.



To prepare the application:

- Cut 13 pads squares;
- Put some skin preparation gel in a cup.
- Put some conductive paste in a cup.
- Put the gloves on.

Applying the IceCap Neonate

To apply the **IceCap Neonate**:

1. Place the electrode Cz on the corresponding electrode position on baby's head. Check the positioning of other electrodes to make sure the IceCap Neonate is placed correctly.

Warning
Do not apply excessive pressure on the electrodes.


Warning
Make sure that the conductive side is in contact with patient's skin. Labels shall be on the outer side.

2. Apply skin preparation gel to the scalp under the Cz electrode.
3. Using your fingers, apply some conductive paste under the Cz electrode.
4. Apply skin preparation gel and conductive paste under the other electrodes.
5. If additional fixation for electrode is needed, apply some adhesive cream on each pad square and place each pad square on top of electrode to help fixing the electrode on the skin.

**Note**

It is possible to use some medical tape and/or a hair net to maintain the electrodes.

6. Plug the **IceCap Neonate** to the IceAdapter or Touchproof Adapter.

**Warning**

Fix the connectors using medical tape or a hair net to avoid any risk of strangulation


**Warning**

Make sure that the L and R marks from the flat cable correspond to the figure above.

**Note**

Other compatible devices are listed in [Section 2.4, “Compatible Devices”](#).

7. Plug the IceAdapter or Touchproof Adapter to a compatible device.
8. Using a compatible device, check the impedance. If the values are too high, remove the hair under the corresponding electrodes and softly push the electrodes against the scalp.

4. SAFETY INFORMATION

4.1. Operating Conditions

The devices are designed to be used in a controlled medical environment (hospital, doctor's practice, private hospital), in a research environment, or at home.

The following environmental conditions must be met in order to use the product safely:

Characteristics	Environmental conditions
Usage temperature	Between +5 °C and +40 °C
Storage temperature	Between +5 °C and +50 °C
Transportation temperature	Between -25 °C and +70 °C
Relative humidity	Between 5% and 95%
Atmospheric pressure	Between 700 hPa and 1060 hPa

4.2. Safety Information

Warning

The IceCaps are electrical devices containing type BF applied parts in contact with the patient in accordance with the standard 60601-1.

Warning

Electromedical devices require special precautions regarding EMC (Electromagnetic Compatibility) and must be installed and operated in accordance with the EMC information provided in the sections below.

Warning

The use of accessories, transducers and cables other than those specified (transducers and cables sold by the manufacturer of the **IceCaps) such as replacement parts for internal components, may cause an increase of emissions or a reduction of immunity of the device.**

Warning

The **IceCaps are intended to be used on still patients.**

**Warning**

Do not use the system simultaneously with any other medical device that has not been validated along with (external defibrillator).

**Warning**

Do not use a damaged or degraded device (packaging damaged, product damaged, corrosion).

**Warning**

Do not modify the device.

**Warning**

Electrodes, conductive parts and associated connectors for applied parts, including the neutral electrode, must not be in contact with other conductive parts, including ground.

4.3. Information on Residual Risks

**Warning**

Do not use the device on a patient with wounds.

**Warning**

It is recommended to use the device away from any extreme source of heat and cigarettes.

**Warning**

Use the system indoor in a quiet environment.

**Warning**

The device will be damaged by water.

**Warning**

To avoid contamination risk of the patient, the healthcare professional shall ensure that the device is clean before using it on a patient.

**Warning**

Do not use the device in case of allergy to one of the components.

**Warning**

Use on healthy skin only.

5. REGULATORY INFORMATION

5.1. System Life Cycle

5.1.1. Product Lifetime

The components of the system lifetime is estimated in the following list:

Component	Number of uses	Shelf life
IceAdapter	200	-
Touchproof Adapter	200	-
IceCap 2	1	2 years
IceCap Neonate	1	2 years

Once the lifetime duration of the product is reached, contact the support.

5.1.2. Recycling the IceCap 2 Parts

Directive 2012/19/EU of the European Parliament and of the Council of 4 July 2012 on waste electrical and electronic equipment (WEEE) and Directive 2006/66/EC of the European Parliament and Council of 6 September 2006 on batteries and accumulators and waste batteries and accumulators and repealing Directive 91/157/EEC governing the recycling of batteries, electrical wire and electronic components.

The electrical wires and other electronic components cannot be discarded along with household waste.

The IceCaps should be discarded and disposed of in an appropriate container.

**Note**

Return the device to the manufacturer for recycling.

5.2. Product Symbology

Symbol	Description
	Manufacturer
	Date of manufacture
	Single use
	Use by date
	Serial number
	Catalogue reference
	Do not use if package is damaged
	Keep away from sunlight
	Lower and upper limit of temperature
	Keep dry
	Upper and lower limit of humidity
	The product must be eliminated in a suitable salvage and recycling structure
	Type BF applied part
	Operating instructions
	Quantity of product
	Size of the product
	Medical device
IP20	IP rating
	Atmospheric pressure limitation
	Refer to instruction manual

Do not expose the device to magnetic resonance environments

Recyclable parts

The device is collected separately

5.2.1. Manufacturer Information

BioSerenity Medical Devices Group
6/8, rue Jean Antoine de Baïf,
75013 Paris - FRANCE

The **IceCap 2**, **IceCap 2 Small** and **IceCap Neonate** are CE marked medical devices, according to regulation (UE) 2017/745.

5.3. Medical Device Classification

EU

The **IceCaps** are class I medical devices.

US and Canada

The **IceCaps** are class II medical devices.

5.4. Contact

5.4.1. Contact Information

To contact the support services, use the following email address: support.neuro@bioserenity.com

5.4.2. Internet

The information provided is available and updated on the website: <https://www.bioserenity.com/>

5.4.3. Instructions for Use

The paper version of the instructions for use is available upon request at no additional cost under 7 calendar days. To receive the instructions for use, contact the support.

5.4.4. Materiovigilance

Every serious incident related to the device must be notified to the manufacturer and to the competent authority from the user's State.

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BioSerenity Medical Device Group - SAS